ZENIT

Light Cure Nano Ceramic Composite

PREBOND

ETCHING GEL

Dynamic*Flow*
Universal Nano Flow Composite

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USER GUIDE

Nano Ceramic Composite is a light-curing, high-lustre-polish hybrid composite with an ultrafine, radiopaque porcelain filler for use in adhesive filling treatment. It can be polished to a high lustre. Due to the ultra-fine particle filler, extremely homogeneous restorations can be placed which are easily polished to a high lustre. The chameleon effect matches the shade of the filling perfectly to the tooth structure. The guidelines in DIN EN ISO 4049 have been complied with. Nano Ceramic Composite is supplied in practical screw syringes.

Composition:

Monomer matrix: diurethane dimethacrylate, butanediol dimethacrylate, isopropylidene-bis[2(3)-hydroxy-3(2)-(4-phenoxy)propyl] bismethacrylate [Bis-GMA]

Total filler content 83% by weight (70% by volume)

- Glass filler (medium grit size 0.7 microns)
- Pyrogenic silica (medium grit size 12 nm)
- Agglomerated nanoparticles (medium grit size 0.6 microns)

Indications

- Direct anterior and posterior restorations in black's class I, II, III, IV and V cavities.
- Indirect restorations such as inlays, onlays and laminate veneers
- Extended fissure sealing in molars and premolars
- Cores
- Splinting mobile teeth
- Adjusting the contours and shade to improve aesthetics

Application

Preparatory Measures

Before commencing the treatment, clean the tooth with non-fluoride polishing paste. Use a Vita shade guide to select the shade while the tooth is still moist.

1. Cavity preparation prepare the cavity minimally invasively as generally required for adhesive techniques All enamel margins in the anterior region must be bevelled. Do not bevel the margins in the posterior region and avoid slice preparations. Spray the cavity with water to clean it, remove all debris and dry it. The cavity must be isolated. It is advisable to place a rubber dam.
2. Pulp protection / Cavity liner. If an enamel-dentine adhesive is used, no cavity liner is required. In very deep cavities those areas in close proximity to the pulp must be coated with a calcium hydroxide material.
3. Approximal contact areas when filling cavities with approximal sections, place a transparent matrix and fix it in place.
4. Etching with Etching Apply Etching first to the enamel areas of the cavity and let it take effect for 30 seconds. Then fill the whole cavity (dentine) with etching gel and let it take effect for a further 15 seconds. The etching time in the dentine should not exceed 20 seconds. Then rinse the cavity thoroughly with a water spray and dry it with oil-free compressed air. Do not overdry the dentine. Dried, etched enamel surfaces have a chalky-white appearance and must not be contaminated before the bonder is applied. If the surface becomes contaminated with saliva, rinse and dry again and re-etch if necessary.
5. Application of the bonding agent refer to the Prebond LC instructions for details.
6. Placing the composite in the cavity remove the required amount of composite from the screw syringe, insert it into the cavity metal instruments and contour. Layers should not be more than 2 mm thick. A dispersion layer, a thin, using standard non-polymerised film, remains on the surface after curing due to the effect of the oxygen in the air. This produces the chemical bond between the layers and must not be touched or contaminated with moisture. Curing the exposure time for all shades per layer is 40 seconds with a standard halogen polymerisation unit or an LED polymerisation lamp or 2x3 seconds with a plasma polymerisation unit. Hold the light guide as close as possible to the surface of the filling. Polymerise each surface of complex fillings.
7. Trimming Nano Ceramic Composite can be trimmed and polished immediately after curing using finishing diamonds, flexible disks, silicone polishers and polishing brushes. Check the occlusion and articulation and spot grind to eliminate high spots or undesirable paths of articulation from the surface of the filling.

Indirect method

Cavity preparation

The cavity should be prepared as minimally invasively as possible with only slightly diverging sides. All internal line and point angles must be rounded. Avoid slice preparations. Prepare a flat cervical shoulder-do not bevel it. Any unavoidable undercuts must be blocked out with glass ionomer cement. Use slightly tapering diamonds with rounded tips for the preparation. Coat those areas of dentine in close proximity to the pulp with a thin layer of calcium hydroxide material. Cavity liners containing eugenol are contraindicated.

Fabricating an inlay

Cast the impression with high strength dental stone in the laboratory. Allow the model to set and pull off the impression. Block out the under-cuts and apply an oil-free separating agent to the model. Build up the inlay layer-by-layer on the model. Build up the approximal and deep occlusal sections first. Each individual layer must not be more than 2 mm thick and is cured separately with a commercially available light curing lamp. The finished inlay is then released from the die and cured fully. Trim the occlusal surface with fissure burs and polish to a high lustre with silicone polishers and diamond paste. Clean the inlay thoroughly with soap and water, spray with air/water to rinse and dry.

Fitting inlays, onlays or veneers

Remove the temporary restoration and clean the cavity. Place a rubberdam, then clean and dry the prepared surface of the tooth. Check the fit of the restoration applying light pressure. Do not force the restoration to fit. Grind the fitting surface if necessary to improve the fit. The occlusion and articulation should not be checked during the try-in of the inlay, as there is a risk of fracturing the inlay. Etch the tooth with Etching (refer to the Etching instructions for use), rinse thoroughly with the water spray and dry the etched surfaces with oil-free compressed air. The dried, etched enamel surface has a chalky-white appearance and must not be contaminated before applying Prebond LC. If the surface is contaminated with saliva, rinse and dry again and re-etch if necessary. Apply a thin layer of Prebond LC bonding agent with a brush to the etched enamel and cavity walls, rub it in and cure for 40 seconds with a standard halogen polymerisation unit. After curing, a dispersion layer forms that should not be removed, as it produces the chemical bond to the filling material. Sandblast the interior of the composite restoration surface and clean with alcohol, then apply Prebond LC without light curing. Mix commercial dual-hardening composite cement and apply to the interior surface of the restoration. Press the inlay carefully into position and remove larger amounts of excess with a spatula. Remove approximal excess with a probe or dental floss. Press lightly on the inlay with a ball-end stopper until curing is complete to ensure that the inlay remains in the correct position.

Beginning at the approximal, cure each surface of the composite for 40 seconds with a standard halogen polymerisation unit or LED polymerisation unit. Remove excess with finegrit diamond rotary instruments and then diamond coated finishing strips. Check the occlusion and adjust it if necessary. Complete final finishing and polishing with a finishing and polishing set.

Note: If an inlay is more than 2 mm thick, a dual-curing composite must be used

Please note

- The working time beneath a dental light is approximately 2 minutes.
- When placing time consuming restorations, to prevent the composite curing prematurely the dental light should be moved away from the site temporarily or the composite covered with foil impervious to light.
- Use a light curing unit with an emission spectrum of 350 – 500 nm for curing this material. As the required physical properties can only be achieved if the lamp is functioning correctly, its luminous intensity must be checked regularly as described by the manufacturer.

Storage

Do not store above 25°C. Avoid direct sunlight. Close the screw syringes tightly immediately after use. The material should be at room temperature before use. Retract the plunger of the syringe slightly to prevent the apertures becoming blocked. Do not use after the expiry date (refer to label on syringe). For use by dentists only. Keep out of reach of children. This product was developed specifically for the described range of applications. It must be used as described in the instructions. The manufacturer is not liable for damage caused by handling or processing the material incorrectly.

Side-effects

With proper use of this medical device, unwanted side-effects are extremely rare. Reactions of the immune system (allergies) or local discomfort, however, cannot be ruled out completely. Should you learn about unwanted side-effects - even if it is doubtful that the side-effect has been caused by our product – please kindly contact us. To prevent possible reactions of the pulp in cavities where the dentine is exposed, the pulp must be protected adequately (e.g. calcium hydroxide preparation).

Contra-indications/interactions

If a patient has known allergens against or hypersensitivities towards a component of this product, we recommend not to use it or to do so only under strict medical supervision. In such cases, we will supply the composition of our medical device upon request. The dentist should consider known interactions and cross reactions of the product with other materials already in the patient's mouth before using the product. Unpolymerized composite may cause skin allergies. The user must take a adequate

precautions. As phenolic substances (such as eugenol) inhibit polymerization, do not use cavity liners (such as zinc-oxide eugenol cements) containing such substances. Vita is a registered trademark of Vita®Zahnfabrik H. Rauter GmbH, & Co.KG, Bad Säckingen, Germany.

Trouble shooting

Problem	Cause	Remedy
Composite does not cure completely	Light output of the light-curing lamp inadequate	Check the light output and change the light source if required
	Emitted wave length range of the light-curing lamp inadequate	Consult the manufacturer of the light-curing lamp. Recommended wavelength range: 350 - 500 nm
Composite in the syringe is sticky and soft, colourless liquid separates in the syringe	Material has been stored over a longer period at $\geq 25^{\circ}\text{C}$	Adhere to storage temperature; store in a refrigerator
	Material has been kept too long in a syringe warmer	Never keep a syringe longer than one hour per application in a syringe warmer
Composite appears too hard and firm in the syringe	The material has not been heated to room temperature after removal from the refrigerator	Allow the composite to heat to room temperature before use; use a syringe warmer if necessary
	Syringe not properly sealed, composite partially cured	Always seal the syringe properly with the cap after removing composite
Inlay/onlay is not properly retained when fitted	Restoration is too opaque to be cemented using only light-curing composite	Use dual-curing luting composite
Composite does not cure all the way through (dark or opaque shades)	Composite layers applied too thickly for each curing cycle	Adhere to a max. layer thickness of 2.0 mm per increment
Restoration appears too yellow compared with the shade guide	Inadequate curing of the composite layer	Repeat the exposure cycle several times; min. 40 sec.

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USER GUIDE

Prebond is an adhesive with outstanding adhesion to dentin and etched enamel. It improves adhesion between the composite and tooth surface.

Indication:

Prebond is used to improve the adhesion of composites on enamel and dentin. Further, it improves the adhesion between glass-ionomer underfill (lining) cement and composite fillings ("sandwich technique").

Application:

The cavity should be prepared according to the rules for adhesive filling technique. Depending on the situation, apply calcium hydroxide paste or an underfill (lining) cement in deep cavities. The enamel is etched, rinsed and then dried in accordance with the etch bonding technique.

The glass-ionomer cement filling should not be etched, only roughened.

Once the prepared cavity is dry, 1-2 drops of light-curing Prebond are dripped onto the mixing pad (immediately reclose the bottle) and a thin layer is applied to the tooth using a brush or foam pellet. Excess material is removed using a dry brush or a pellet. The material is then cured using halogen light (15 sec.).

Side effects:

Undesirable side-effects are extremely rare if the product is prepared and used properly. Should you become aware of any undesirable side-effects, including in cases of doubt, please notify us, giving as precise a description of the attendant circumstances and

symptoms as possible.

Working instructions:

Ensure that the working instructions and precautionary measures described in the supplied product information are always strictly observed. Should the patient suffer undesirable reactions or any other unwanted incidents arise even though these instructions were observed, then please inform us and give a precise description of the accompanying circumstances.

Contra-indications:

Do not use Prebond in the area of the pulpa.

Storage/shelf-life:

3 years, store in a cool, dry place. Protect against light!

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USER GUIDE

Description

Etching Gel is phosphoric acid formulation thickened with natural polymeric materials. Etching Gel contains 37% phosphoric acid and is designed as a semi-gel type. Etching Gel can be obtained easily and available etching effect. Also Etching Gel leaves no residues on the surface of etched tooth contrary to silica thickened etchants.

Composition

- Phosphoric acid
- Purified water
- Thickener
- Colorant

Indications

Etching the enamel and dentine for adhesive restorations.

Contra-indications

Do not apply on soft tissues including oral mucosa, eyes and dental pulp tissue.

Recommended Procedure

1. Isolate tooth and prepare the cavity in a conservative manner.
2. Equip the disposable tip after removing the cap.
3. Etch enamel and dentin with etching agent for 15 seconds. Rinse gel thoroughly and remove excess water from the preparation with a gentle stream of air for 1-2 seconds.

4. Dry etched surfaces with clean air. Do not blow saliva onto etched surfaces. If etched surface does not appear chalky white, repeat etching. Salivary contamination must be avoided.
5. Repeat steps 3-4 if desired.

Storage

1. Close the syringe cap after use.
2. Do not store at elevated temperature or intense light.

Precautions

1. For dental use only.
2. It should be used by a dentist and dental professionals. Observe the procedure and the precautions.
3. Keep it away from reach of non-dental users including children, senior citizens and others.
4. Avoid contacting with eye and skin
5. If contacts skin, wash immediately with soap and water.
6. Using protective gloves is recommended.
7. If accidental contact with eyes or prolonged oral tissues occurs, flush immediately with a large amount of water and consult a physician.
8. If ingestion, give a large amount of water or milk, take an antacid, seek medical attention.
9. Do not use after expiry date.

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Farbenton: A1, A2, A3, A3.5, A4, B1, B2, B3, Oa2, OA3, C2, D2, WT, TL

USER GUIDE

Dynamic Flow is a low viscosity, light-cured, radiopaque flowable composite resin for invasive cavity preparation.

Indication

Restoration of invasive cavity

Pit and fissure sealing

Base and Liner for direct restoration

Repair acrylic temporary materials

Repair of esthetic indirect restorations

Contra-Indications

Dynamic Flow should not be used if patients are known to be allergic to methacrylate resins.

Instructions

Before isolating and preparing choose the desired shade and type of Dynamic Flow. Clean with a slurry of pumice and use and dry.

Apply 37% PD etching gel to available enamel dentin 15 sec., remove gel and rinse for 20 sec. Blot dry, keep moist and do not rub.

Apply bond to the tooth. Apply double coat of President Prebond by using microbrush. Air dry gently for 5 sec. Light Cure for 10

sec. if you apply other products, use them according to manufacturers instructions. composite Dispensing-Dynamic Flow flowable

can be delivered directly from the disposable tip. Place and light cure restorative in increments. Cure each increment by its entire surface to curing light for 20 seconds. Finish and Polish-After contouring composite, polish surface with finishing and polishing proximal areas. Occlusion-Check lateral and centric occlusion adjust as necessary.

Cautions:

- Avoid contact of Dynamic Flow with the skin/mucous membrane or eyes.
- None polymerized Dynamic Flow may causes lightly, reversible irritation or lead to a general sensitization against methacrylates.
- Keep away from children. For dental use only.
- Do not use after the date of expiration.
- Exposure to light may lead to premature polymerization.

Storage

Do not store above 22 °C After use close syringe cap.

Avoid from direct sun light.

Shelf Life 3 years after manufacturing date

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PRESIDENT DENTAL GmbH