DE

# ZIRCON.X MULTILAYER®

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### Instruction For Use

#### Before Use

President Dental Multilayer zirconia blocks are produced by CIP technology and pre-sintered in low temperature. Although the product has some strength, because of porosity, please handle carefully. When you receive the product, please check as below. If there is something exceptions, please contact President Dental GmbH.

- Product is complete without any damage
- The packing is complete without any damage
- Label: Company name, product name, batch number, inspectors and inspection date.

#### [Multilayer-3D Material Usage and Technical Parameters] Material Function

The zirconia embryo used in dentistry, can be used to manufacture crown, the bridges and other restoration.

#### Chemical Composition and Powder Characteristics

##### Chemical Composition

Zr <sub>2</sub> R <sub>3</sub>	90%-95%
Y <sub>2</sub> O <sub>3</sub>	4%-10%
Al <sub>2</sub> O <sub>3</sub>	≤0.5%
Other Oxide	<0.5%
Aging Properties	Monoclinic phase <5%
Chemical Solubility	≤2000µg.cm <sup>-1</sup>
Sintered Chemical Solubility	≤100µg.cm <sup>-1</sup>

##### Mechanical Properties

Density	≥2.8 g/cm <sup>3</sup>
Sintered Density	≥6.0g/cm <sup>3</sup>
Flexural Strength	Inzisal 600MPa Cervical 800MPa
Fracture Toughness	≥3.5Mpam <sup>-1/2</sup>
Sintered Hardness (Hv10)	Hv10±12.5Gpa

#### [Application Range]

Veneer, Inlays, Onlays, Anterior Crown, Posterior Crown, Full Contour Bridge (Below 6 Units)

#### [Color]

A1,A2,A3,A3.5,A4 B1, B2, B3, B4 C1, C2, C3, C4 D2, D3, D4

#### [Requirements for Preparation]


Multilayer zirconia is not suitable for non-shoulder with slant. Since both ways will lead to a thinner edge and risk of fracture.

#### [Anterior Crown:

The inner edge of shoulder must be slick or be fluted.

Preparation width of incisal, labial and palatal should be above 0.7mm.

The axial wall must be blunt to ensure the zirconia can reach an ideal milling effect with CAD/CAM.


#### [Full Contour:

The inner edge of shoulder must be slick or be fluted. The minimum width of incisal edge, labial side, and palatal side is at least 0.7mm.

Width of occlusal surface of lateral side of the lip, tongue preparation is more than 1.0mm

The axial wall must be blunt to ensure the zirconia can reach an ideal milling effect with CAD/CAM.

#### [Application Method]

Full contour crown and bridges:

Scanning → Designing → Milling → Cleaning → Sintering → Polishing → Staining → Finishing

#### 1. Scanning

**Oral Scanning:** Before oral scanning, the first thing is to do gingival treatment to expose the cervical line and clean and dry the blood and saliva in scanning area.

The initial scanning region is the preparation region, from the maxilla region to the buccal lingual side, and the the jaw.  
Note: when scanning the edge, the speed should be slow and repeat scanning.

#### Scanning the butment:

Scanning the base teeth one by one. Making sure that the base teeth should make a full contact with model.

#### 2. Designing

**Coping:** The thinnest should be no less than 0.7mm Increase compensation data of the milling Burs. It is an effective way to easy the condition that the anterior teeth cannot match the veneer of butment well due to the keenness of the veneer of butment. The tip of the butment should be filled with wax, and then processing the scanning and designing.


#### 3. Milling

When making restorations by using President Dental zirconia ceramic materials, make sure that a new bur should be used. And to cool zirconia by using liquid is avoided when doing milling process. After milling is finished, check that if there are any defects happening listed below:

Is there any crack? Is there any contamination? Is there any break? If any of these defects happen, reasons must be found and restorations need to be milled again.

#### [Layout:

Choose zirconia block of suitable thickness. When layout, put data in the middle of disc (vertical direction). The distance between tooth data and disc surface should be at least 1mm. Connector should be placed in the most protruding position on buccal side. The thickness of connector should be 2mm. (as picture shows)

#### [Place the disc:

Directing arrow is seen on the side of "Multilayer" disc. The arrow points at the direction of incisal, namely the part of disc showing shallower color.

#### 1. Separation of Restoration and Cleaning

##### Separation of restoration:

Using technician specialized hand piece and grinding head to separate restorations from blocks.

Before grinding, a towel should be put on the desk to avoid restorations dropping on desk and crack or break. When operating, hand needs to find a fulcrum: hand piece speed should be controlled at 10000-12000 rev/min. And then the connector need to be polished successively in one direction. Dont separated one connector completely off at one time; Finally, the rest of connector could be polished slightly. It is not advised to make too much adjustment on the restorations in soft condition to avoid causing subsfissure or chipping and so on.

##### Cleaning:

Use President Dental special zirconia grinding head to trim the surface of restoration. The following three procedures (coarse grinding, fine grinding and coarse polishing) can make restoration surface smooth.

**Coarse grinding:** This is the first grinding step after sintering restoration. The aim is to seat the restoration, adjust adjoining, occlusion and trim the anatomic contour of teeth.

**Fine grinding:** Please maket he tooth surface evenly and until 7, ormlly as well as the surface texture more smoothly.

#### 4. Staining and Glazing

Please stain the entire ceramic crown or entire zirconia crown.

#### 5. Finished

Fabrication of restoration is finished.

#### Attachment

The suggested crystallization curve (10h):

Start Temp	Room Temp
Heating time in the first stage	90 min
Max. Temp. in the first stage	1100°C
Heating time in the second stage	250min
Max.Temp. in the second stage	1400°C
Holding Time	10min
Heating time in the third stage	50min
Max. Temp in the third stage	1500°C
Holding time	120min
Cooling time	70min
Cooling to	900°C

**Stain & glaze sintering process:** (To avoid subsfissure on bridge, the heating rate for bridge stain & glaze sitering must go down to 30°C/min)

Start Temp	500°C
Dry	4mins
Rate	50°C/min
Highest Temp	820°C
Maintain	2mins
Final Temp.	300°C

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